



CQ-03140-D2O-COSY-271009-298 K-500 MHz  
 COSY-mmcd-clem D2O /opt/topspin ravi {1 A4 - 104}



NAME CQ-03140-D2O-271009  
 EXPNO 8  
 PROCNO 1  
 Date\_ 20091029  
 Time 10.53  
 INSTRUM spect  
 PROBHD 5 mm CPTXO 13C  
 PULPROG cosyqf45  
 TD 2048  
 SOLVENT D2O  
 NS 8  
 DS 16  
 SWH 6493.506 Hz  
 FIDRES 3.170657 Hz  
 AQ 0.1577460 sec  
 RG 11.3  
 DW 77.000 usec  
 DE 30.00 usec  
 TE 298.0 K  
 D0 0.00000300 sec  
 D1 2.00000000 sec  
 IN0 0.00015400 sec

ppm

-1

0

1

2

3

4

5

6

7

8

9

10

11

===== CHANNEL f1 =====  
 NUC1 1H  
 P1 15.70 usec  
 PL1 0.90 dB  
 PL1W 21.64242172 W  
 SFO1 499.8423492 MHz  
 ND0 1  
 TD 256  
 SFO1 499.8423 MHz  
 FIDRES 25.365242 Hz  
 SW 12.991 ppm  
 FnMODE QF  
 SI 1024  
 SF 499.8400000 MHz  
 WDW SINE  
 SSB 0  
 LB 0.00 Hz  
 GB 0  
 PC 1.40  
 SI 1024  
 MC2 QF  
 SF 499.8400000 MHz  
 WDW SINE  
 SSB 0  
 LB 0.00 Hz  
 GB 0

ppm